

5 questions for your healthcare provider

Learn whether a smart¹ insulin pen is for you

1. How might a smart¹ insulin pen help improve my life with diabetes?

2. Who can access my dosing information?

3. How can my dosing information be used during visits with my healthcare provider and in my daily life?

4. How do I get the most out of tracking my diabetes dosing with a smart¹ insulin pen?

5. How do I get started with a smart¹ insulin pen?



What makes an insulin pen smart¹?

NovoPen® 6 is a smart¹ insulin pen. What makes it smart¹ is that it records your insulin dosing information automatically, so you don't have to note it down yourself. And the pen allows you to easily check information about your last dose, so you don't have to worry about remembering it.

You can transfer this information to your preferred compatible app*, allowing you to view your dosing information alongside glucose information from your Blood Glucose Monitor or Continuous Glucose Monitor (if you are using one). Plus, all of this information can be shared with your healthcare provider or family members, if you wish.

NovoPen® 6))

*Compatible apps vary by region. Speak with your healthcare provider to find out what is currently available for your NovoPen® 6.

For people living with diabetes who require the administration of insulin.

Reference: 1. Sangave NA, et al. *Diabetes Spectr* 2019;32:378–84.

Novo Nordisk Pharmaceuticals Pty Ltd. Level 10, 118 Mount Street, North Sydney NSW 2060. ABN 40 002 879 996. Novo Nordisk® Medical Information (Australia) 1800 668 626. www.novonordisk.com.au. ® Registered trademark of Novo Nordisk A/S. NONOY0060. AU24DI00042. October 2024.